


Packing Instructions

STP-309DI Biological Substance, Category B Insulated Shipping System



Packing Instructions

1. Follow the manufacturer's closure instructions to close the primary receptacle(s). Ensure that all primary receptacle closures are held securely, tightly and effectively in place by secondary means. For example, tape may be used around the closure.
2. Complete address for shipper and consignee.
3. Assemble the outer packaging and secure the bottom with packing tape. Place insulated chest inside the outer packaging.
4. For multiple fragile primary receptacles, use the provided STP-609 Saf-T-Rap to separate and protect the primary receptacles. Place the primary receptacle(s) in the leakproof polybag portion (STP-701) of the patented Secondary Pressure Receptacle (STP-700) along with absorbent material. Follow the directions on the polybag to seal the bag.

NOTE: Use an amount of absorbent sufficient to absorb the entire volume of the primary receptacle(s). The absorbent included in the system will absorb up to 50 mL.

5. Complete an Itemized List of Contents directly on the STP-700 Tyvek® Envelope with an indelible pen.
6. Place the sealed polybag into the Tyvek® Envelope and seal according to the instructions on the envelope.
7. Place the completed Secondary Pressure Receptacle into the STP-309 inner box. Use void fill for space within the inner box. Place the closed inner box into the bottom of the insulated chest.

NOTE: For primary receptacles with a volume greater than 50 mL, the closures of the primary receptacle(s) must be oriented with the orientation arrows on the outer packaging.

8. Place up to 2 kg of dry ice in the insulated chest, around and on top of the inner box.
9. Place lid onto the insulated chest. Close the fibreboard flaps and secure with packing tape.

For performance data with Dry ice, please contact Saf-T-Pak at:
1-800-814-7484 from the US or Canada or 780-486-0211 worldwide.

Capacity and Payload Options



Test Tubes

Maximum Capacity:

Eight (8), 10 mL tubes.
Four tubes per STP-700.
One STP-700 provided.

NOTE: The maximum payload capacity of the STP-309DI is 95 g or 95 mL.
The maximum weight capacity of the STP-309DI is 3 kg.